

WNT2 Antibody FITC Conjugated

Catalog No: #C06238F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	WNT2 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 240-290 360 derived from human WNT2
Conjugates	FITC
Target Name	WNT2
Other Names	IRP; INT1L1; Protein Wnt-2; Int-1-like protein 1; Int-1-related protein; WNT2
Accession No.	Swiss-Prot#P09544NCBI Gene ID7472
Uniprot	P09544
GeneID	7472;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters.

Note: This product is for in vitro research use only