

RAC1+RAC2 Antibody HRP Conjugated

Catalog No: #C06248H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RAC1+RAC2 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 115-145 192 derived from human RAC1 RAC2.
Conjugates	HRP
Target Name	RAC1+RAC2
Other Names	MIG5; p21 Rac1; Rac 1; Rac-1; p21 Rac2; PAK1; Ras like protein TC25; Ras related C3 botulinum toxin substrate 1; Ras related C3 botulinum toxin substrate 2; RAC1_HUMAN; RAC2_HUMAN.
Accession No.	NCBI Gene ID5879, 5880
GeneID	58795880
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

The protein encoded by this gene is a GTPase which belongs to the RAS superfamily of small GTP-binding proteins. Members of this superfamily appear to regulate a diverse array of cellular events, including the control of cell growth, cytoskeletal reorganization, and the activation of protein kinases. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2009].

Note: This product is for in vitro research use only