

FOXP3 Monoclonal Antibody

Catalog No: #27179

Package Size: #27179-1 50ul #27179-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	FOXP3 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	7H9-D6-A10
Isotype	IgG1
Applications	WB ICC
Species Reactivity	Hu
Specificity	Transfected Only
Immunogen Type	Recombinant Protein
Immunogen Description	Purified recombinant human FOXP3 protein fragments expressed in E.coli.
Target Name	FOXP3
Other Names	AIID; AIID; DIETER; DIETER; Forkhead box P3; Forkhead box protein P3; foxp3; foxp3; FOXP3_HUMAN; FOXP3delta7; Immune dysregulation polyendocrinopathy enteropathy X linked; Immunodeficiency polyendocrinopathy enteropathy X nked; IPEX; IPEX; JM2; MGC141961
Accession No.	Uniprot: Q9BZS1 Gene ID: 50943
Uniprot	Q9BZS1
GeneID	50943;
SDS-PAGE MW	47kd
Formulation	Ascites containing with 0.02% sodium azide and 50% glycerol.
Storage	store at -20Λ C

Application Details


Western blotting: 1:300

Immunocytochemistry: 1:100

Images



Immunocytochemistry staining of CHO-K1 cells which transfected by FOXP3 fixed in 1% Paraformaldehyde and using anti-FOXP3 antibody (dilution 1:100).



Western blot detection of FOXP3 antibody in CHO-K1 cell lysate. A0 and CHO-K1 transfected by FOXP3. B0 cell lysate using FOXP3 antibody (1:300 diluted). Predicted band size: 47KDa. Observed band size: 50KDa.

Background

Probable transcription factor. Plays a critical role in the control of immune response.

Note: This product is for in vitro research use only