

CYP2W1 Antibody HRP Conjugated

Catalog No: #C06326H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CYP2W1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human CYP2W1
Conjugates	HRP
Target Name	CYP2W1
Other Names	CP2W1_HUMAN; CYP2W 1; CYP2W1; CYP11W1; Cytochrome P450 2W1; Cytochrome P450, family 2, subfamily W, polypeptide 1; MGC34287.
Accession No.	NCBI Gene ID54905
Uniprot	Q8TAV3
GeneID	54905;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids.

Note: This product is for in vitro research use only