

## FUT8 Antibody FITC Conjugated

Catalog No: #C06329F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FUT8 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 320-370 575 derived from human FUT8
Conjugates	FITC
Target Name	FUT8
Other Names	Alpha-(1,6)-fucosyltransferase; Alpha1-6FucT; Fucosyltransferase 8; GDP-L-Fuc:N-acetyl-beta-D-glucosaminide alpha1,6-fucosyltransferase; GDP-fucose--glycoprotein fucosyltransferase; Glycoprotein 6-alpha-L-fucosyltransferase; FUT8
Accession No.	Swiss-Prot#Q9BYC5NCBI Gene ID2530
Uniprot	Q9BYC5
GeneID	2530;
Excitation Emission	494nm 518nm
Cell Localization	Intracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

## Background

Catalyzes the addition of fucose in alpha 1-6 linkage to the first GlcNAc residue, next to the peptide chains in N-glycans.

Note: This product is for in vitro research use only