

FUT8 Antibody HRP Conjugated

Catalog No: #C06329H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FUT8 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 320-370 575 derived from human FUT8
Conjugates	HRP
Target Name	FUT8
Other Names	Alpha-(1,6)-fucosyltransferase; Alpha1-6FucT; Fucosyltransferase 8; GDP-L-Fuc:N-acetyl-beta-D-glucosaminide alpha1,6-fucosyltransferase; GDP-fucose--glycoprotein fucosyltransferase; Glycoprotein 6-alpha-L-fucosyltransferase; FUT8
Accession No.	Swiss-Prot#Q9BYC5NCBI Gene ID2530
Uniprot	Q9BYC5
GeneID	2530;
Excitation Emission	N A
Cell Localization	Intracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

Catalyzes the addition of fucose in alpha 1-6 linkage to the first GlcNAc residue, next to the peptide chains in N-glycans.

Note: This product is for in vitro research use only