

## TMPRSS4 Antibody FITC Conjugated

Catalog No: #C06334F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TMPRSS4 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 80-130 437 derived from human TMPRSS4
Conjugates	FITC
Target Name	TMPRSS4
Other Names	CAPH2; MT-SP2; TMPRSS3; Transmembrane protease serine 4; Channel-activating protease 2; Membrane-type serine protease 2; TMPRSS4; UNQ776 PRO1570
Accession No.	Swiss-Prot#Q9NRS4NCBI Gene ID56649
Uniprot	Q9NRS4
GeneID	56649;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

## Background

Probable protease. Seems to be capable of activating ENaC (By similarity).

Note: This product is for in vitro research use only