

## Annexin VI Monoclonal Antibody

Catalog No: #27198

Package Size: #27198-1 50ul #27198-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Annexin VI Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	8A7-F8-G12
Isotype	IgG1
Purification	Affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	This antibody detects endogenous levels of ANXA6 and does not cross-react with related proteins.
Immunogen Type	Recombinant Protein
Immunogen Description	Purified recombinant human ANXA6 protein fragments expressed in E.coli.
Target Name	Annexin VI
Other Names	67 kDa calelectrin; Annexin VI; Annexin-6; Calphobindin-II; Chromobindin-20; Lipocortin VI; Protein III; p68; p70; ANXA6;
Accession No.	Uniprot: P08133 Gene ID: 309
Uniprot	P08133
GeneID	309;
SDS-PAGE MW	76kd
Formulation	Purified mouse monoclonal in PBS(pH 7.4) containing with 0.2% sodium azide, 50% glycerol, 0.1%BSA
Storage	store at -20Λ C

## Application Details

Western blotting: 1:1000

## Images



Western blot detection of ANXA6 in T47D, MB4, CEM, MOLT-4 & Jurkat cell lysates and using ANXA6 antibody (1:1000 diluted). Predicted band size: 76KDa Observed band size: 76KDa

## Background

---

May associate with CD21. May regulate the release of Ca<sup>2+</sup> from intracellular stores.

---

Note: This product is for in vitro research use only