

## SMYD4 Monoclonal Antibody

Catalog No: #27215

Package Size: #27215-1 50ul #27215-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	SMYD4 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	7C8-B11-F10
Isotype	IgG1
Applications	WB
Species Reactivity	Hu
Specificity	This antibody detects endogenous levels of SMYD4 and does not cross-react with related proteins.
Immunogen Type	Recombinant Protein
Immunogen Description	Purified recombinant human SMYD4 protein fragments expressed in E.coli.
Target Name	SMYD4
Other Names	KIAA1936; SET and MYND domain containing 4; SET and MYND domain-containing protein 4; SMYD4; SMYD4_HUMAN; ZMYND21.
Accession No.	Uniprot: Q8IYR2 Gene ID: 114826
Uniprot	Q8IYR2
GeneID	114826;
SDS-PAGE MW	89kd
Formulation	ascites
Storage	store at -20 $\wedge$ C


## Application Details

Western blotting: 1:200

## Images



Western blot detection of SMYD4 antibody in 40ug HeLa, 40ug MCF7 and 40ug MDA-MB-468 cell lysates using SMYD4 antibody (1:100 diluted). Predicted band size: 89kDa. Observed band size: 89kDa.



Western blot detection of SMYD4 antibody in HeLa, MCF7 and A431 cell lysates using SMYD4 antibody (1:200 diluted). Predicted band size: 89KDa. Observed band size: 89KDa.

Note: This product is for in vitro research use only