

Dlx1 Antibody HRP Conjugated

Catalog No: #C06428H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dlx1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human Dlx1
Conjugates	HRP
Target Name	Dlx1
Other Names	DII B; Distal less homeo box 1; Distal less homeobox 1; Distal-less homeobox gene 1a; Dlx 1; Dlx; Dlx1; Dlx1 distal-less homeobox 1; DLX1_HUMAN; Homeo box protein DLX1; Homeobox protein DLX 1; Homeobox protein DLX-1; Homeobox protein Dlx1a; Homeodomain transcription factor DLX1; RP23-419H3.3.
Accession No.	NCBI Gene ID1745
Uniprot	P56177
GeneID	1745;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Likely to play a regulatory role in the development of the ventral forebrain. May play a role in craniofacial patterning and morphogenesis and may be involved in the early development of diencephalic subdivisions.

Note: This product is for in vitro research use only