

## CKAP4 Antibody FITC Conjugated

Catalog No: #C06449F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CKAP4 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 552-602 602 derived from human CKAP4
Conjugates	FITC
Target Name	CKAP4
Other Names	p63; CLIMP-63; ERGIC-63; Cytoskeleton-associated protein 4; 63-kDa cytoskeleton-linking membrane protein; CKAP4
Accession No.	Swiss-Prot#Q07065NCBI Gene ID10970
Uniprot	Q07065
GeneID	10970;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

## Background

High-affinity epithelial cell surface receptor for APF. Mediates the anchoring of the endoplasmic reticulum to microtubules.

Note: This product is for in vitro research use only