

SKALP Antibody FITC Conjugated

Catalog No: #C06457F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SKALP Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human SKALP Elafin
Conjugates	FITC
Target Name	SKALP
Other Names	ELAF_HUMAN; Elafin; Elastase Specic Inhibitor; Elastase-specic inhibitor; ES 1; ESI; Peptidase inhibitor 3; PI 3; PI-3; PI3; Protease inhibitor 3 skin derived; Protease inhibitor WAP3; SKALP; Skin derived Anti leukoproteinase; Skin-derived antileukoproteinase; Trappin 2; WAP four disulfide core doma
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

SKALP

Note: This product is for in vitro research use only