

BMP15 Antibody HRP Conjugated

Catalog No: #C06507H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | BMP15 Antibody HRP Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | WB |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic peptide derived from human BMP15 |
| Conjugates | HRP |
| Target Name | BMP15 |
| Other Names | BMP 15; BMP-15; BMP15; BMP15_HUMAN; Bone morphogenetic protein 15; GDF 9B; GDF-9B; GDF9B; Growth differentiation factor 9B; ODG2; POF4. |
| Accession No. | NCBI Gene ID9210 |
| Uniprot | O95972 |
| GeneID | 9210; |
| Excitation Emission | N A |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:500-2000

Background

The BMP15 protein is a member of the bone morphogenetic protein family which is part of the transforming growth factor beta superfamily. The transforming growth factor beta superfamily includes large families of growth and differentiation factors. It is thought that BMP15 may be involved in oocyte maturation and follicular development as a homodimer, or by forming heterodimers with a related protein, Gdf9.

Note: This product is for in vitro research use only