

Neurotensin Antibody HRP Conjugated

Catalog No: #C06517H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurotensin Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 14-161 170 derived from human Neurotensin(C-QLYENKPRRPYIL)
Conjugates	HRP
Target Name	Neurotensin
Other Names	NN; NT; NT N; NTS1; NMN-125; Neurotensin neuromedin N; NTS
Accession No.	Swiss-Prot#P30990NCBI Gene ID4922
Uniprot	P30990
GeneID	4922;
Excitation Emission	N A
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

Neurotensin may play an endocrine or paracrine role in the regulation of fat metabolism. It causes contraction of smooth muscle.

Note: This product is for in vitro research use only