

IL-10 Antibody HRP Conjugated

Catalog No: #C06590H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-10 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 127-159 178 derived from human IL-10
Conjugates	HRP
Target Name	IL-10
Other Names	CSIF; TGIF; GVHDS; IL-10; IL10A; Interleukin-10; Cytokine synthesis inhibitory factor; IL10
Accession No.	Swiss-Prot#P22301NCBI Gene ID3586
Uniprot	P22301
GeneID	3586;
Excitation Emission	N A
Cell Localization	Secreted, Helical
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Inhibits the synthesis of a number of cytokines, including IFN-gamma, IL-2, IL-3, TNF and GM-CSF produced by activated macrophages and by helper T-cells.

Note: This product is for in vitro research use only