

MSR1 Antibody HRP Conjugated

Catalog No: #C06591H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MSR1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 230-272 451 derived from human MSR1 CD204
Conjugates	HRP
Target Name	MSR1
Other Names	SRA; SR-A; CD24; phSR1; phSR2; SCARA1; Macrophage scavenger receptor types I and II; Macrophage acetylated LDL receptor I and II; Scavenger receptor class A member 1; MSR1
Accession No.	Swiss-Prot#P21757NCBI Gene ID4481
Uniprot	P21757
GeneID	4481;
Excitation Emission	N A
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Membrane glycoproteins implicated in the pathologic deposition of cholesterol in arterial walls during atherogenesis. Two types of receptor subunits exist. These receptors mediate the endocytosis of a diverse group of macromolecules, including modified low density lipoproteins (LDL). Isoform III does not internalize acetylated LDL.

Note: This product is for in vitro research use only