

Aquaporin 8 Antibody HRP Conjugated

Catalog No: #C06604H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aquaporin 8 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 129-156 261 derived from human Aquaporin 8
Conjugates	HRP
Target Name	Aquaporin 8
Other Names	AQP-8; Aquaporin-8; AQP8
Accession No.	Swiss-Prot#O94778NCBI Gene ID343
Uniprot	O94778
GeneID	343;
Excitation Emission	N A
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

Forms a water-specific channel; mercury-sensitive.

Note: This product is for in vitro research use only