

SLC26A4 Antibody FITC Conjugated

Catalog No: #C06605F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLC26A4 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 317-344 780 derived from human Solute carrier family 26 member 4
Conjugates	FITC
Target Name	SLC26A4
Other Names	EVA; PDS; DFNB4; TDH2B; Pendrin; Sodium-independent chloride iodide transporter; Solute carrier family 26 member 4; SLC26A4
Accession No.	Swiss-Prot#O43511NCBI Gene ID5172
Uniprot	O43511
GeneID	5172;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Background

Sodium-independent transporter of chloride and iodide.

Note: This product is for in vitro research use only