

## 14-3-3 zeta+delta Antibody PE Conjugated

Catalog No: #C06607P

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	14-3-3 zeta+delta Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 51-100 245 derived from human 14-3-3 zeta delta
Conjugates	PE
Target Name	14-3-3 zeta+delta
Other Names	HEL4; YWHAD; KCIP-1; HEL-S-3; 14-3-3-zeta; 14-3-3 protein zeta delta; Protein kinase C inhibitor protein 1; YWHAZ
Accession No.	Swiss-Prot#P63104NCBI Gene ID7534
Uniprot	P63104
GeneID	7534;
Excitation Emission	480,565nm 578nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

## Background

Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner.

Note: This product is for in vitro research use only