

## MAGEC1 Antibody FITC Conjugated

Catalog No: #C06631F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MAGEC1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide derived from human MAGEC1
Conjugates	FITC
Target Name	MAGEC1
Other Names	Cancer testis antigen 7.1; Cancer testis antigen family 7 member 1; CT7; CT7.1; MAGC1_HUMAN; MAGE C1; MAGE C1 antigen; MAGE-C1 antigen; MAGEC1; melanoma antigen family C, 1; Melanoma associated antigen C1; Melanoma-associated antigen C1; MGC39366.
Accession No.	NCBI Gene ID9947
Uniprot	O60732
GeneID	9947;
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

## Background

MAGEC1

Note: This product is for in vitro research use only