

## HFREP1 Antibody HRP Conjugated

Catalog No: #C06672H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HFREP1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human HFREP1 FGL1
Conjugates	HRP
Target Name	HFREP1
Other Names	FGL 1; Fgl1; FGL1_HUMAN; Fibrinogen like 1; Fibrinogen like protein 1; Fibrinogen related protein 1; Fibrinogen-like protein 1; Hepassocin; Hepatocellular carcinoma related sequence; Hepatocyte derived fibrinogen related protein 1; Hepatocyte-derived fibrinogen-related protein 1; HFREP 1; HFREP-1; H
Accession No.	NCBI Gene ID2267
Uniprot	Q08830
GeneID	2267;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200

## Background

Fibrinogen-like 1 is a member of the fibrinogen family. This protein is homologous to the carboxy terminus of the fibrinogen beta- and gamma- subunits which contains the four conserved cysteines of fibrinogens and fibrinogen related proteins. However, this protein lacks the platelet-binding site, cross-linking region and a thrombin-sensitive site which are necessary for fibrin clot formation. This protein may play a role in the development of hepatocellular carcinomas. Four alternatively spliced transcript variants encoding the same protein exist for this gene. [provided by RefSeq, Jul 2008].

Note: This product is for in vitro research use only