

CD81 Antibody FITC Conjugated

Catalog No: #C06702F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD81 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Dog
Immunogen Description	KLH conjugated synthetic peptide aa 112-162 236 derived from human TAPA1 CD81
Conjugates	FITC
Target Name	CD81
Other Names	S5.7; CVID6; TAPA1; TSPAN28; CD81 antigen; 26 kDa cell surface protein TAPA-1; Target of the antiproliferative antibody 1; Tetraspanin-28; Tspan-28; CD81
Accession No.	Swiss-Prot#P60033NCBI Gene ID975
Uniprot	P60033
GeneID	975;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test IF=1:50-200

Background

May play an important role in the regulation of lymphoma cell growth. Interacts with a 16-kDa Leu-13 protein to form a complex possibly involved in signal transduction. May act as the viral receptor for HCV.

Note: This product is for in vitro research use only