

Neutrophil Elastase Antibody PE Conjugated

Catalog No: #C06723P

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neutrophil Elastase Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 100-140 267 derived from human Neutrophil Elastase
Conjugates	PE
Target Name	Neutrophil Elastase
Other Names	GE; NE; HLE; HNE; ELA2; SCN1; PMN-E; Neutrophil elastase; Bone marrow serine protease; Elastase-2; Human leukocyte elastase; Medullasin; PMN elastase; ELANE
Accession No.	Swiss-Prot#P08246NCBI Gene ID1991
Uniprot	P08246
GeneID	1991;
Excitation Emission	480,565nm 578nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Background

Modifies the functions of natural killer cells, monocytes and granulocytes. Inhibits C5a-dependent neutrophil enzyme release and chemotaxis.

Note: This product is for in vitro research use only