

Neutrophil Elastase Antibody HRP Conjugated

Catalog No: #C06723H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Neutrophil Elastase Antibody HRP Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | Flow-Cyt |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic peptide aa 100-140 267 derived from human Neutrophil Elastase |
| Conjugates | HRP |
| Target Name | Neutrophil Elastase |
| Other Names | GE; NE; HLE; HNE; ELA2; SCN1; PMN-E; Neutrophil elastase; Bone marrow serine protease; Elastase-2; Human leukocyte elastase; Medullasin; PMN elastase; ELANE |
| Accession No. | Swiss-Prot#P08246NCBI Gene ID1991 |
| Uniprot | P08246 |
| GeneID | 1991; |
| Excitation Emission | N A |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

Flow-Cyt=1:50-200

Background

Modifies the functions of natural killer cells, monocytes and granulocytes. Inhibits C5a-dependent neutrophil enzyme release and chemotaxis.

Note: This product is for in vitro research use only