

PRAME Antibody FITC Conjugated

Catalog No: #C06734F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PRAME Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms
Immunogen Description	KLH conjugated synthetic peptide aa 180-227 509 derived from human PRAME
Conjugates	FITC
Target Name	PRAME
Other Names	MAPE; OIP4; CT130; OIP-4; Melanoma antigen preferentially expressed in tumors; Opa-interacting protein 4; Preferentially expressed antigen of melanoma; PRAME
Accession No.	Swiss-Prot#P78395NCBI Gene ID23532
Uniprot	P78395
GeneID	23532;
Excitation Emission	494nm 518nm
Cell Localization	Nucleus, Cell membrane
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Background

Functions as a transcriptional repressor, inhibiting the signaling of retinoic acid through the retinoic acid receptors RARA, RARB and RARG. Prevents retinoic acid-induced cell proliferation arrest, differentiation and apoptosis.

Note: This product is for in vitro research use only