

## Interferon alpha 11 Antibody HRP Conjugated

Catalog No: #C06744H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interferon alpha 11 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human IFNA11 Interferon alpha 11
Conjugates	HRP
Target Name	Interferon alpha 11
Other Names	Limitin; N alpha 11; N-alpha 11; na11; NAB_MOUSE; Interferon alpha-11; N-alpha-11; interferon alpha-1.
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

## Background

Produced by macrophages, IFN alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.

Note: This product is for in vitro research use only