

NUPR1 Antibody HRP Conjugated

Catalog No: #C06801H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NUPR1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human NUPR1 p8
Conjugates	HRP
Target Name	NUPR1
Other Names	Candidate of metastasis 1; COM1; Nuclear protein 1; NUPR1; NUPR1_HUMAN; p8 protein; Protein p8.
Accession No.	NCBI Gene ID26471
Uniprot	O60356
GeneID	26471;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

Could participate in the response to proapoptotic stimuli and promotes cellular growth in a way that helps the tissue counteract diverse injuries. May contribute to the metastatic phenotype. Tissue specificity: Highly expressed in pancreas.

Note: This product is for in vitro research use only