

SPON1 Antibody FITC Conjugated

Catalog No: #C06859F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SPON1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 775-807 807 derived from human SPON1 F spondin
Conjugates	FITC
Target Name	SPON1
Other Names	f-spondin; VSGP F-spondin; Spondin-1; Vascular smooth muscle cell growth-promoting factor; SPON1; KIAA0762; VSGP
Accession No.	Swiss-Prot#Q9HCB6NCBI Gene ID10418
Uniprot	Q9HCB6
GeneID	10418;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Cell adhesion protein that promotes the attachment of spinal cord and sensory neuron cells and the outgrowth of neurites in vitro. May contribute to the growth and guidance of axons in both the spinal cord and the PNS (By similarity). Major factor for vascular smooth muscle cell.

Note: This product is for in vitro research use only