

RHBDD3 Antibody FITC Conjugated

Catalog No: #C06891F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RHBDD3 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human RHBDD3 (101-200aa)
Conjugates	FITC
Target Name	RHBDD3
Other Names	C22orf3;HS984G1A;Pituitary tumor apoptosis;PTAG;RHBD3_HUMAN;RHBDD3;Rhomboid domain containing 3;Rhomboid domain-containing protein 3.
Accession No.	NCBI Gene ID25807
Uniprot	Q9Y3P4
GeneID	25807;
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

RHBDD3

Note: This product is for in vitro research use only