

RHCG Antibody FITC Conjugated

Catalog No: #C06924F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RHCG Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 30-71 479 derived from human RHCG
Conjugates	FITC
Target Name	RHCG
Other Names	RHGK; PDRC2; C15orf6; SLC42A3; Ammonium transporter Rh type C; Rh glycoprotein kidney; Rhesus blood group family type C glycoprotein; Rh family type C glycoprotein; Rh type C glycoprotein; Tumor-related protein DRC2; RHCG; CDRC2
Accession No.	Swiss-Prot#Q9UBD6NCBI Gene ID51458
Uniprot	Q9UBD6
GeneID	51458;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular, Helical
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Functions as an electroneutral and bidirectional ammonium transporter. May regulate transepithelial ammonia secretion.

Note: This product is for in vitro research use only