

GPR65 Antibody FITC Conjugated

Catalog No: #C06939F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPR65 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 72-89 337 derived from human GPR65
Conjugates	FITC
Target Name	GPR65
Other Names	TDAG8; hTDAG8; Psychosine receptor; G-protein coupled receptor 65; T-cell death-associated gene 8 protein; GPR65
Accession No.	Swiss-Prot#Q8IYL9NCBI Gene ID8477
Uniprot	Q8IYL9
GeneID	8477;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Receptor for the glycosphingolipid psychosine (PSY) and several related glycosphingolipids. May have a role in activation-induced cell death or differentiation of T-cells.

Note: This product is for in vitro research use only