

Anillin Antibody FITC Conjugated

Catalog No: #C06984F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anillin Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 503-548 1124 derived from human Anillin
Conjugates	FITC
Target Name	Anillin
Other Names	scra; FSGS8; Scraps; Actin-binding protein anillin; ANLN
Accession No.	Swiss-Prot#Q9NQW6NCBI Gene ID54443
Uniprot	Q9NQW6
GeneID	54443;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm, Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Background

Required for cytokinesis. Essential for the structural integrity of the cleavage furrow and for completion of cleavage furrow ingression.

Note: This product is for in vitro research use only