

## Factor X heavy chain Antibody FITC Conjugated

Catalog No: #C07003F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Factor X heavy chain Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human Factor X heavy chain
Conjugates	FITC
Target Name	Factor X heavy chain
Other Names	Activated factor Xa heavy chain; Coagulation factor; Coagulation factor X; EC 3.4.21.6; F10 antibody FA10_HUMAN; Factor Xa; FX; FXA; Prothrombinase; Stuart factor; Stuart Prower factor; Stuart-Prower factor.
Accession No.	NCBI Gene ID2159
Uniprot	P00742
GeneID	2159;
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

## Background

Factor Xa is a vitamin K-dependent glycoprotein that converts prothrombin to thrombin in the presence of factor Va, calcium and phospholipid during blood clotting.

Note: This product is for in vitro research use only