

NUBP1 Antibody HRP Conjugated

Catalog No: #C07042H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NUBP1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human NUBP1
Conjugates	HRP
Target Name	NUBP1
Other Names	Cytosolic Fe S cluster assembly factor NUBP1; Cytosolic Fe-S cluster assembly factor NUBP1; NBP 1; NBP; NBP1; nubp1; NUBP1_HUMAN; Nucleotide binding protein e.coli MinD like; Nucleotide binding protein 1 E.coli MinD like; Nucleotide binding protein 1 MinD homolog, E. coli; Nucleotide binding protein
Accession No.	NCBI Gene ID4682
Uniprot	P53384
GeneID	4682;
Excitation Emission	N A
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Implicated in the regulation of centrosome duplication (By similarity). Component of the cytosolic iron-sulfur (Fe S) protein assembly machinery. Required for maturation of extramitochondrial Fe S proteins. May bind and transfer 2 labile 4Fe-4S clusters to target apoproteins.

Note: This product is for in vitro research use only