

RABGAP1 Antibody FITC Conjugated

Catalog No: #C07047F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RABGAP1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human RABGAP1
Conjugates	FITC
Target Name	RABGAP1
Other Names	GAP and centrosome associated protein; HSPC094; Rab GTPase activating protein 1; Rab6 GTPase activating protein GAPCenA; RBGP1_HUMAN.
Accession No.	NCBI Gene ID23637
Uniprot	Q9Y3P9
GeneID	23637;
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

May act as a GTPase-activating protein of RAB6A. May play a role in microtubule nucleation by centrosome. May participate in a RAB6A-mediated pathway involved in the metaphase-anaphase transition.

Note: This product is for in vitro research use only