

Cyclin T2 Antibody HRP Conjugated

Catalog No: #C07089H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin T2 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human Cyclin T2
Conjugates	HRP
Target Name	Cyclin T2
Other Names	CCNT2; CCNT2_HUMAN; Cyclin-T2; CycT2; FLJ90560; MGC134840.
Accession No.	NCBI Gene ID905
Uniprot	O60583
GeneID	905;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Regulatory subunit of the cyclin-dependent kinase pair (CDK9 cyclin T) complex, also called positive transcription elongation factor B (P-TEFB), which is proposed to facilitate the transition from abortive to production elongation by phosphorylating the CTD (carboxy-terminal domain) of the large subunit of RNA polymerase II (RNAP II).

Note: This product is for in vitro research use only