

GAB2 (Phospho-Ser159) Antibody

Catalog No: #11771

Package Size: #11771-1 50ul #11771-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GAB2 (Phospho-Ser159) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of GAB2 only when phosphorylated at serine 159.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around phosphorylation site of Serine 159(K-S-S(p)-A-P) derived from Human GAB2.
Target Name	GAB2
Modification	Phospho
Other Names	GABJ; KIAA0571; pp100;
Accession No.	Swiss-Prot#: Q9UQC2; NCBI Gene#: 9846; NCBI Protein#: NP_536739.1.
Uniprot	Q9UQC2
GeneID	9846;
SDS-PAGE MW	100kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500~1:1000

Images



Western blot analysis of extracts from HT-29 cells treated with serum using GAB2 (Phospho-Ser159) Antibody #11771. The lane on the right is treated with the antigen-specific peptide.

Background

This gene is a member of the GRB2-associated binding protein (GAB) gene family and is similar to the GAB1 gene. These proteins contain pleckstrin homology (PH) domain, and bind SHP2 tyrosine phosphatase and GRB2 adapter protein. They act as adapters for transmitting various signals in response to stimuli through cytokine and growth factor receptors, and T- and B-cell antigen receptors. The protein encoded by this gene is the principal activator of phosphatidylinositol-3 kinase in response to activation of the high affinity IgE receptor. Two alternatively spliced transcripts encoding different isoforms have been described for this gene.

Nishida K., *Blood* 93:1809-1816(1999).

Nagase T., *DNA Res.* 5:31-39(1998).

Taylor T.D., *Nature* 440:497-500(2006).

Note: This product is for in vitro research use only