

PID1 Antibody FITC Conjugated

Catalog No: #C07200F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PID1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms
Immunogen Description	KLH conjugated synthetic peptide derived from human PID1
Conjugates	FITC
Target Name	PID1
Other Names	P CLI1; Ab1 119; P-CLI1; PCLI1; PCLI1_HUMAN; Phosphotyrosine interaction domain-containing protein 1; PID1; Protein NYGGF4; PTB containing, cubilin and LRP1 interacting protein; PTB-containing, cubilin and LRP1-interacting protein.
Accession No.	NCBI Gene ID55022
Uniprot	Q7Z2X4
GeneID	55022;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Increases proliferation of preadipocytes without affecting adipocytic differentiation. Tissue specificity: Expressed in subcutaneous fat, heart, skeletal muscle, brain, colon, thymus, spleen, kidney, liver, small intestine, placenta, lung and peripheral blood leukocyte.

Note: This product is for in vitro research use only