

Alpha casein Antibody FITC Conjugated

Catalog No: #C07336F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha casein Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide derived from human Alpha casein
Conjugates	FITC
Target Name	Alpha casein
Other Names	Alpha casein; CASA; CASA1_HUMAN; Casein alpha s1; Casoxin-D; CSN1; CSN1S1.
Accession No.	NCBI Gene ID1446
Uniprot	P47710
GeneID	1446;
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Important role in the capacity of milk to transport calcium phosphate. Casoxin D acts as opioid antagonist and has vasorelaxing activity mediated by bradykinin B1.

Note: This product is for in vitro research use only