

DPY19L2 Antibody FITC Conjugated

Catalog No: #C07371F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DPY19L2 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human DPY19L2
Conjugates	FITC
Target Name	DPY19L2
Other Names	D19L2_HUMAN; dpy 19 like 2; dpy 19 like 2 C. elegans; Dpy 19 like protein 2; Dpy-19-like protein 2; dpy19; DPY19L2; FLJ32949; FLJ36166; Protein dpy 19 homolog 2; Protein dpy-19 homolog 2.
Accession No.	NCBI Gene ID283417
Uniprot	Q6NUT2
GeneID	283417;
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Required during spermatogenesis for sperm head elongation and acrosome formation.

Note: This product is for in vitro research use only