

BEND3 Antibody HRP Conjugated

Catalog No: #C07482H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BEND3 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 530-575 828 derived from human BEND3
Conjugates	HRP
Target Name	BEND3
Other Names	KIAA1553; BEN domain-containing protein 3; BEND3
Accession No.	Swiss-Prot#Q5T5X7NCBI Gene ID57673
Uniprot	Q5T5X7
GeneID	57673;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

BEND3

Note: This product is for in vitro research use only