

## FNDC5 Antibody PE Conjugated

Catalog No: #C07514P

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FNDC5 Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human FNDC5 Irisin
Conjugates	PE
Target Name	FNDC5
Other Names	FRCP2; irisin; Fibronectin type III domain-containing protein 5; Fibronectin type III repeat-containing protein 2; FNDC5
Accession No.	Swiss-Prot#Q8NAU1NCBI Gene ID252995
Uniprot	Q8NAU1
GeneID	252995;
Excitation Emission	480,565nm 578nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

## Background

Irisin: Contrary to mouse, may not be involved in the beneficial effects of muscular exercise, nor in the induction of browning of human white adipose tissue.

Note: This product is for in vitro research use only