

Pancreatic Polypeptide Antibody HRP Conjugated

Catalog No: #C07550H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pancreatic Polypeptide Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 30-65 95 derived from human Pancreatic hormone
Conjugates	HRP
Target Name	Pancreatic Polypeptide
Other Names	PP; PNP; Pancreatic prohormone; Pancreatic polypeptide; PPY
Accession No.	Swiss-Prot#P01298NCBI Gene ID5539
Uniprot	P01298
GeneID	5539;
Excitation Emission	N A
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Pancreatic hormone is synthesized in pancreatic islets of Langerhans and acts as a regulator of pancreatic and gastrointestinal functions. The physiological role for the icosapeptide has not yet been elucidated.

Note: This product is for in vitro research use only