

TCL1A Antibody

Catalog No: #32079

Package Size: #32079-1 50ul #32079-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TCL1A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total TCL1A protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human TCL1A (NP_001092195.1).
Target Name	TCL1A
Other Names	TCL1A;TCL1
Accession No.	Uniprot:P56279GeneID:8115
Uniprot	P56279
GeneID	8115
SDS-PAGE MW	16kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of Ramos cells, using TCL1A antibody.

Background

Overexpression of the TCL1 gene in humans has been implicated in the development of mature T cell leukemia, in which chromosomal rearrangements bring the TCL1 gene in close proximity to the T-cell antigen receptor (TCR)-alpha (MIM 186880) or TCR-beta (MIM 186930) regulatory elements (summarized by Virgilio et al., 1998 [PubMed 9520462]). In normal T cells TCL1 is expressed in CD4⁺/CD8⁻ cells, but not in cells at later stages of differentiation. TCL1 functions as a coactivator of the cell survival kinase AKT (MIM 164730) (Laine et al., 2000 [PubMed 10983986]).[supplied by OMIM, Jul 2010]

Note: This product is for in vitro research use only