

TCTN3 Antibody

Catalog No: #32081

Package Size: #32081-1 50ul #32081-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TCTN3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total TCTN3 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human TCTN3 (NP_056446.4).
Target Name	TCTN3
Other Names	TCTN3;C10orf61;JBTS18;OFD4;TECT3
Accession No.	Uniprot:Q6NUS6GeneID:26123
Uniprot	Q6NUS6
GeneID	26123
SDS-PAGE MW	68kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TCTN3 antibody.

Background

This gene encodes a member of the tectonic gene family which functions in Hedgehog signal transduction and development of the neural tube. Mutations in this gene have been associated with Orofaciodigital Syndrome IV and Joubert Syndrom 18. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

Note: This product is for in vitro research use only