

ATG13 Antibody

Catalog No: #32084

Package Size: #32084-1 50ul #32084-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ATG13 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total ATG13 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human ATG13 (NP_055556.2).
Target Name	ATG13
Other Names	KIAA0652;PARATARG8;ATG13
Accession No.	Uniprot:O75143GeneID:9776
Uniprot	O75143
GeneID	9776
SDS-PAGE MW	75kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ATG13 antibody.

Background

The protein encoded by this gene is an autophagy factor and a target of the TOR kinase signaling pathway. The encoded protein is essential for autophagosome formation and mitophagy.

Note: This product is for in vitro research use only