

## Somatostatin Antibody HRP Conjugated

Catalog No: #C07641H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Somatostatin Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 75-107 116 derived from human Somatostatin
Conjugates	HRP
Target Name	Somatostatin
Other Names	SMST; Somatostatin; Growth hormone release-inhibiting factor; SST
Accession No.	Swiss-Prot#P61278NCBI Gene ID6750
Uniprot	P61278
GeneID	6750;
Excitation Emission	N A
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200

## Background

Somatostatin inhibits the release of somatotropin.

Note: This product is for in vitro research use only