

AKIRIN2 Antibody HRP Conjugated

Catalog No: #C07712H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AKIRIN2 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human AKIRIN2
Conjugates	HRP
Target Name	AKIRIN2
Other Names	AKIR2_HUMAN; AKIRIN 2; Akirin-2; Akirin2; C6orf166; dJ486L4.2; FBI1; FLJ10342; Fourteen three three beta interactant 1.
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

Required for the innate immune response. Downstream effector of the Toll-like receptor (TLR), TNF and IL-1 beta signaling pathways leading to the production of IL-6. Forms a complex with YWHAB that acts to repress transcription of DUSP1.

Note: This product is for in vitro research use only