

FBXO27 Antibody HRP Conjugated

Catalog No: #C07716H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FBXO27 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human FBXO27
Conjugates	HRP
Target Name	FBXO27
Other Names	F-box G-domain protein 5; F-box only protein 27; F-box protein 27; F-box protein FBG5; FBG5; Fbx27.
Accession No.	NCBI Gene ID126433
Uniprot	Q8NI29
GeneID	126433;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Members of the F-box protein family, such as FBXO27, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1, cullin (CUL1), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains.

Note: This product is for in vitro research use only