

TYMP Antibody

Catalog No: #32154

Package Size: #32154-1 50ul #32154-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TYMP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total TYMP protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human TYMP (NP_001107228.1).
Target Name	TYMP
Other Names	TYMP; ECGF; ECGF1; MEDPS1; MNGIE; MTDPS1; PDECGF; TP; hPD-ECGF
Accession No.	Uniprot:P19971 GeneID:1890
Uniprot	P19971
GeneID	1890
SDS-PAGE MW	50kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TYMP antibody.

Background

This gene encodes an angiogenic factor which promotes angiogenesis in vivo and stimulates the in vitro growth of a variety of endothelial cells. It has a highly restricted target cell specificity acting only on endothelial cells. Mutations in this gene have been associated with mitochondrial neurogastrointestinal encephalomyopathy. Multiple alternatively spliced transcript variants have been identified.

Note: This product is for in vitro research use only